

DECLARATION OF COMPLIANCE Seitz® T Series Filter Sheets


Product Grades :

T120, T500, T750, T950, T1000, T1500, T2100, T2600, T3500, and T5500

The Seitz T series filter sheets are comprised of a cellulose matrix with perlite. They also contain resin and polyolefin fibers.

An initial flush is recommended prior to use.

Issued	1 March 2010
Revised	29 September 2017
Expires	28 February 2019
Reference	FBDCTENg
Page	1 of 2



Dawn Cromwell
Director, QA
Laboratory, Food & Beverage

Seitz T Series Filter Sheets

Components

Cellulose, binder resin, perlite and polyolefin fibers

Declaration

Europe

Seitz T series filter sheets meet the requirements of food contact as detailed in European Regulation (EC) Number 1935/2004 in that:

- Our supplier states that the monomers and additives of the polymer fibers are listed in European Regulation (EU) Number 10/2011 Annex I.
- The product complies with German Recommendations XXXVI and XXXVI/1, as well as with the German Foodstuffs and Animal Feed Code (LFGB) §§ 30 and 31.

USA

The following materials of construction meet the FDA requirements for food contact use as detailed in the US Code of Federal Regulations, 21 CFR paragraphs 177-199 in that:

- Cellulose and binder resin to 21 CFR section 177.2260 (Filters, resin bonded) and to 21 CFR section 176.170 (Components of paper and paperboard in contact with aqueous and fatty foods)
- Polyolefin fibers to 21 CFR section 177.1520 (Olefin polymers)
- Total extractables as per 21 CFR 177.2260 (Filters, Resin Bonded) (g) (h) (i) (j) (k) (l).
50% ethanol at room temperature and n-hexane at reflux were used in the extractables testing.

The following are listed in the Food Chemical Codex (FCC): Perlite

Process Quality System

Site of Manufacture: Pall Filtersystems GmbH, Bad Kreuznach, Germany on behalf of Pall International Sàrl.

The Quality Management System at Pall Filtersystems GmbH, Bad Kreuznach is certified to ISO 9001:2008. These products / product packaging carry a lot number / date code to facilitate traceability to suppliers' materials and Pall production records.

Pall Filtersystems GmbH confirm that the product manufacturing environment, for the above product at our Bad Kreuznach site, is in line with the principles of food contact materials GMP as detailed in Regulation 2023/2006.

Supplied in Europe by

Pall International Sàrl
Av. de Tivoli 3
Fribourg
Switzerland
CH-1700




Pall Food and Beverage

New York - USA
+1 516 484 3600 telephone
+1 866 905 7255 toll free
foodandbeverage@pall.com

Visit us on the web at www.pall.com

Pall Corporation has offices throughout the world. For Pall representatives in your area, please go to www.pall.com/contact.

Because of developments related to products, systems and/or services described herein, the data and procedures are subject to change without notice. Please consult your Pall representative or visit www.pall.com to verify that this information remains valid.

© Copyright 2017, Pall Corporation. Pall,  and Seitz are trademarks of Pall Corporation. © indicates a trademark registered in the USA. Filtration.Separation.Solution is a service mark of Pall Corporation.

FBDCTENg

September 2017