



STATEMENT OF COMPLIANCE: ISO 12500-3:2009 Emflon® PFRW Filter Cartridges AB Style “W” Code

Cartridge Part Number



This is a guide to the part numbering structure only. For specific options, please contact Pall.

Table 1 : Nominal Length

Code	Description
05*	127 mm (5")
1	254 mm (10")
2	508 mm (20")
3	762 mm (30")

*Available in codes 2 and 7 only

Table 2 : Adaptor

Code	Description
2	SOE – single open end with flat closed end, 2 locking tabs and external 226 O-rings
3	SOE – single open end with flat closed end and external 222 O-rings
7	SOE – single open end with fin end, 2 locking tabs and external 226 O-rings
8	SOE – single open end with fin end and external 222 O-rings
28	SOE – single open end with fin end, 3 locking tabs and external 222 O-rings

Table 3 : Seal Material

Code	Description
H4	Silicone elastomer
J	Ethylene propylene rubber

Emflon PFRW final air filters are designed for gas filtration in the food and beverage industry, fermenter inlet air and exhaust venting, sterile process air and sterile venting of tanks in food and beverage processing.

'Sterile air' is defined as that filtered using a 0.2 µm liquid rated filter (capable of removing 10⁷ of *Brevundimonas diminuta* per cm² of effective filtration area, in laboratory tests). Please note this product is not supplied pre-sterilized.

The key feature of this product is the inclusion of two layers of hydrophobic, 0.2 µm microbial rated Emflon PTFE medium in a polypropylene hardware filter construction.

Issued 26 February 2013
 Revised 11 April 2017
 Expires 28 February 2019
 Page 1 of 2

Steven Bailey
 Quality Manager
 Pall Manufacturing (UK) Ltd., Ilfracombe

Emflon PFRW Filter Cartridges (AB Style “W” Code)

Components

Filter Media	Two layers of hydrophobic, Emflon PTFE media
Support / Drainage	Polypropylene
Side Seal	Polypropylene
Cage, Core, End Cap and Fin End	Polypropylene
Adaptor	Polypropylene (with integral stainless steel reinforcing ring as appropriate to design*)
O-ring Seal	Silicone Elastomer (H4) Ethylene Propylene Rubber (J)

*For applications requiring autoclaving or steam sterilization, Pall recommends the use of Code 2 or Code 7 adaptors to ensure filter sealing after cooling. Cartridges should be cooled to system operating temperature prior to use.

Declaration

Emflon PFRW filters meet the requirements of ISO12500-3:2009 in that:

Emflon PFRW filter cartridges (tested as part number: AB1PFR7WH4) have demonstrated applicable performance capability in laboratory tests as described in report LrF 03b/12 appended. These studies were conducted on randomly sampled product made with the same Pall proprietary hydrophobic PTFE filter membrane, polypropylene support/drainage material and polypropylene hardware; and using the same manufacturing techniques.

Supplied in Europe by

Pall International Sàrl
Av. de Tivoli 3
Fribourg
Switzerland
CH-1700




Pall Food and Beverage

New York - USA
+1 516 484 3600 telephone
+1 866 905 7255 toll free
foodandbeverage@pall.com

Visit us on the web at www.pall.com

Pall Corporation has offices throughout the world. For Pall representatives in your area, please go to www.pall.com/contact.

Because of developments related to products, systems and/or services described herein, the data and procedures are subject to change without notice. Please consult your Pall representative or visit www.pall.com to verify that this information remains valid.

© Copyright 2017, Pall Corporation. Pall,  and Emflon are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. Filtration.Separation.Solution is a service mark of Pall Corporation.

April 2017